


**AUTOCLAVE**

**YOUR OFFICE SAFER**

**21**  
SISTEMAS DE SEGURANÇA



**AUTOCLAVE**  
**AUTOCLAVE 21L**

**TECHNOLOGY EVOLVES YOU TOO**

Designed to deliver high performance in all your procedures, only the Saevo Autoclave line integrates technology with design. Guaranteed health and safety at its highest level.

- Intelligent temperature, pressure and altitude control system;
- Cross temperature and pressure data for greater accuracy
- Electronic and mechanical devices that offer greater protection
- Drying cycle with door ajar or door closed
- Bivolt (automatic selection)
- Altitude calibration up to 4,000m for more accurate pressure control

- 01 VANTI-VACUUM VALVE**  
Low pressure
- 02 SAFETY VALVE**  
Overpressure
- 03 SYSTEM OF DEPRESSURIZATION**

- 04 SILICONE SEALING RING**  
Excellent flexibility and resistance to temperature and pressure ranges.

- 05 3 TRAYS**

- 06 STAINLESS STEEL TANK**  
It allows easy and efficient cleaning.

- 07 MANOMETER**  
Allows monitoring of proper temperature and operating pressure throughout the sterilization cycle.

- 08 COMMAND PANEL**  
Command panel located on the side of the equipment. Protecting against temperature changes

**AUTOCLAVE**  
**AUTOCLAVE 12L**



2 BANDEJAS



ANVISA: 10069210082



**COMMAND PANEL**

- 01 | Indication of pressure and temperature
- 02 | Beginning of the sterilization cycle
- 03 | Drying
- 04 | Cancels/Stops the sterilization cycle
- 05 | Equipment connected
- 06 | Indication of heating/sterilization process
- 07 | Indication of the drying process
- 08 | Interrupting cycle

ANVISA - Sealer and Distiller: Non-correlated product. APR - 2024

REQUEST YOUR QUOTE



AUTOCLAVE

imagem ilustrativa.