

DENTAL UNIT

S 200 F

SAEVO

START INSPIRED

The S200 will start your career with technology, quality and design



Illustrative image.

GET TO KNOW S 200

- 1 Biarticulated headrest
- 2 Upholstery in 15 colour options
- 3 Foot control to activate and program the chair
- 4 Anti-skid edging
- 5 No need for floor mounting
- 6 Allows coupling of up to 5 terminals
- 7 Stainless steel tray
- 8 Delivery unit with fixed arm
- 9 Cup-holder and proximity sensor*
- 10 Dental light: Sirius 3, 5 or 3x2 LEDs with proximity sensor*
- 11 Assistant arm with control pad*

VIEW THE LIST OF STANDARD AND OPTIONAL ITEMS ON THE BACK

*Optional in the models S200 F, S200 SF and S200 Cart

DENTAL UNIT S 200 SF



Standard Dental Light: Persus LED

DENTAL UNIT S 200 CART



ASSISTANT ARM WITH CONTROL PAD

UPGRADE: EXTENDED ARM

For greater comfort and safety.



Possibility of including up to three terminals (2 suction units and 1 triple syringe):

- 1x Venturi suction terminal
- 1x VAC PLUS suction terminal*
- 1x External suction terminal*
- 1x Triple Syringe terminal*

PAD with elevation and descent control of the seat and backrest*

*Optional in the models S200 F, S200 SF and S200 Cart

SAEVO COLOURS

Upholstery that inspires your work routine

PVC



LEATHER*



*Optional in the models S200 F, S200 SF and S200 Cart

STANDARD/OPTIONAL

See what comes with your **S 200**

N: NOT APPLICABLE | O: OPTIONAL | S: STANDARD ITEM

DENTAL UNIT	F	SF	CART
Trendelenburg Position	N	N	N
Massage System	N	N	N
PVC Upholstery	S	S	S
Viscoelastic Upholstery	N	N	N
Leather Upholstery	O	O	O
Removable Headrest	N	N	N
Biarticulated Headrest	S	S	S
Armrests: 2 Fixed	S	S	S
Armrests: 1 Fixed + 1 Folding Lateral	O	O	O
Armrests: 1 Fixed + 1 Folding and Pneumatic	N	N	N
Foot Protector	O	O	O
Control Pedal (7 keys)	S	S	S
Joystick Control Pedal (5 keys)	N	N	N
Joystick Pedal (5 keys) - with Chip Blower Function	N	N	N
Foot Control Pedal - with Chip Blower function	N	N	N
Soft Comfort	O	O	O
Connection Box	O	O	O
Ambidextrous	O	O	O

DELIVERY UNIT	F	SF	CART
Pneumatic Lock Arm	S	N	N
Electrical Cables Kit	O	O	O
Stainless Steel Cover	S	S	S
Upper Tray (Acrylic)	N	N	N
Lower Tray (Injected)	N	N	N
Kit PAD (Electronic Controls)**	O	O	O
X-ray Viewer Kit**	O	O	O
Flush Kit	O	O	O
Manometer	N	N	N
High Speed Terminal	S	S	S
Low Speed Terminal - without spray	S	S	S
Low Speed Terminal - with Spray	N	N	N
Triple Syringe	S	S	S
Injected Triple Syringe	N	N	N
Triple Syringe Heating Kit**	O	O	N
Number of terminals in the delivery unit	UP TO 5 (3S & 2O)	UP TO 5 (3S & 2O)	UP TO 5 (3S & 2O)
Electric Micromotor Kit**	O	O	O
Borden Terminal Kit - with Spray	O	O	O
Fiber-Optic Terminal Kit**	O	O	O
Curing Light Kit**	O	O	O
Ultrasonic Scaler Kit**	O	O	O
Bicarbonate Jet Kit	O	O	O

WATER UNIT	F	SF	CART
Cup Holder Kit - with Sensor**	O	O	O
Cup Holder Kit - without Sensor**	O	O	O
Electric Water Bowl Activator	N	N	N
Assistant Arm - without PAD	O	O	O
Assistant Arm - with PAD**	O	O	O
Venturi Suction Kit	S	S	S
Vac Plus Suction Kit	O	O	O
External Suction Terminal Kit	O	O	O
Triple Syringe Kit	O	O	O
Injected Triple Syringe Kit	N	N	N
Triple Syringe Heating Kit**	O	O	O

DENTAL LIGHT	F	SF	CART
Persus LED	S	S	S
Sirius Sensor - 3 LEDs	O	O	O
Sirius Sensor - 5 LEDs	O	O	O
Sirius Sensor - 5 (3x2) LEDs	O	O	O

OTHERS	F	SF	CART
Multimedia Kit Support	O	O	O

**Purchasing the electrical cables kit additionally, is necessary.

DENTAL SET S200 F, S200 SF AND S200 CART CONFIGURATION:

- Dental Chair: Syncrus • Anvisa: 10069210062
- Water Unit: Syncrus G2 • Anvisa: 10069210063
- Delivery Unit: Syncrus G2 • Anvisa: 10069210075
- Dental Light: Persus LED • Anvisa: 10069210069
- Illustrative images • May 2023

REQUEST YOUR QUOTE
S200