

DENTAL UNIT

S 300 F

SAEVO

CHANGE YOUR ROUTINE WITH NEW POSSIBILITIES

The S300 is set to meet your professional demands - with features that will transform your work routine.



Illustrative image.

GET TO KNOW S 300

- 1 Biarticulated headrest
- 2 90° folding arm support
- 3 Upholstery with 15 colour options
- 4 Joystick control pedal
- 5 Anti-skid edging
- 6 No need for floor mounting
- 7 Allows coupling of up to 5 terminals
- 8 Support tray for instruments in version H and stainless steel tray in version F
- 9 Retractable rods with lock (version H)
- 10 Delivery unit with pneumatic arm
- 11 90° folding water unit, cup holder and proximity sensor
- 12 Dental light: Sirius 5 or 3x2 LEDs with proximity sensor*
- 13 Assistant arm with control pad*
- 14 Ergonomic stool*

*Optional.

VIEW THE LIST OF STANDARD AND OPTIONAL ITEMS ON THE BACK

DENTAL UNIT

S 300 H



WATER UNIT



EXTENDED ARM WITH PAD*
*Optional item

Proximity Sensor

Water unit and vat retractable spittoon at 90°

Removable

Water Flow Control of the spittoon and cup filler

Extended Arm* comes with:

- 1x Venturi Suction terminal
- 1x VAC PLUS Suction terminal*
- 1x External Suction terminal*
- 1x Triple Syringe terminal*

*optional

SAEVO COLOURS

Upholstery that inspires your work routine

PVC OR VISCOELASTIC*



LEATHER*



*Optional. Illustrative images

STANDARD/OPTIONAL

See what comes with your **S 300**

N: NOT APPLICABLE | O: OPTIONAL | S: STANDARD ITEM

DENTAL UNIT	F	H
Trendelenburg Position	N	N
Massage System	N	N
PVC Upholstery	S	S
Viscoelastic Upholstery	O	O
Leather Upholstery	O	O
Removable Headrest	N	N
Biarticulated Headrest	S	S
Armrests: 2 Fixed	N	N
Armrests: 1 Fixed + 1 Folding Lateral	S	S
Armrests: 1 Fixed + 1 Folding and Pneumatic	N	N
Foot Protector	S	S
Control Pedal (7 keys)	N	N
Joystick Control Pedal (5 keys)	S	S
Joystick Pedal (5 keys) - with Chip Blower Function	N	N
Foot Control Pedal - with Chip Blower function	N	N
Soft Comfort	O	O
Connection Box	O	O
Ambidextrous	O	O

DELIVERY UNIT	F	H
Pneumatic Lock Arm	S	S
Electrical Cables Kit	S	S
Stainless Steel Cover	S	N
Upper Tray (Acrylic)	N	N
Lower Tray (Injected)	N	S
Kit PAD (Electronic Controls)	O	O
X-ray Viewer Kit	O	O
Flush Kit	O	O
Manometer	N	O
High Speed Terminal	S	S
Low Speed Terminal - without spray	N	N
Low Speed Terminal - with Spray	S	S
Triple Syringe	S	S
Injected Triple Syringe	N	N
Triple Syringe Heating Kit	O	O
Number of terminals in the delivery unit	UP TO 5 (3S & 2O)	UP TO 5 (3S & 2O)
Electric Micromotor Kit	O	O
Borden Terminal Kit - with Spray	O	O
Fibre Optic Terminal Kit	O	O
Curing Light Kit	O	O
Ultrasonic Scaler Kit	O	O
Bicarbonate Jet Kit	O	O

WATER UNIT	F	H
Cup Holder Kit - with Sensor	S	S
Cup Holder Kit - without Sensor	O	O
Electric Water Bowl Activator	S	S
Assistant Arm - without PAD	O	O
Assistant Arm - with PAD	O	O
Venturi Suction Kit	S	S
Vac Plus Suction Kit	O	O
External Suction Terminal Kit	O	O
Triple Syringe Kit	O	O
Injected Triple Syringe Kit	N	N
Triple Syringe Heating Kit	O	O

DENTAL LIGHT	F	H
Persus LED	O	O
Sirius Sensor - 3 LEDs	S	S
Sirius Sensor - 5 LEDs	O	O
Sirius Sensor - 5 (3x2) LEDs	O	O

OTHERS	F	H
Multimedia Kit Support	O	O

DENTAL SET S300 F AND S300 H CONFIGURATION:
 Dental Chair: Syncrus G2 • Anvisa: 10069210062
 Water Unit: Syncrus G2 • Anvisa: 10069210063
 Delivery Unit: Syncrus G2 / Syncrus G3 H • Anvisa: 10069210075
 Dental Light: Sirius G8 Sensor - 3 LEDs • Anvisa: 10069210069
 COD: 51000000210 • May 2023

REQUEST YOUR QUOTE
S300