



ERGONOMICS FOR DEMANDING PROFESSIONALS

Illustrative image.

HANDPIECES

TURBINE HIGH SPEED

SAEVO SL 30 PB ALUMINIUM

- Body made of aluminium
- Rotation of 335.00rpm/min
- Ceramic bearing
- Weight: 30 g
- Low operating noise
- Drill fixation: Push-button
- Torque: 0,05-0,18 N.cm
- * Drill not included

SAEVO SL 30 FG ALUMINIUM

- Body made of aluminium
- Borden connection
- Triple spray
- Rotation: 335.00rpm/min
- Ceramic bearing
- Weight: 30 g
- Low operating noise
- Drill Fixation: Friction Grip
- Torque: 0,05-0,18 N.cm



Fixation of drills
PUSH BUTTON



Triple spray



Fixation of drills
FRICTION GRIP



Borden
and Midwest



Bit breaker



CONTRA-ANGLE

SAEVO SL 30 PB | SAEVO SL 30 FG

- Made of aluminium
- Intra-rotary connection system
- Weight: 49g
- Transmission: 1:1
- Torque: 0,35 to 1,00 N.cm
- * Drill not included.



Fixation of drills
PUSH BUTTON



External Spray

MICROMOTOR

SAEVO SL 30

- Made of aluminium
- Rotation: 5.000 - 20.000 rpm
- Intra-rotary connection system
- Weight: 80g
- Torque: 0,35 to 1,00 N.cm
- Low working noise
- Left and Right direction rotation
- Spray adjustment



With external spray

STRAIGHT PIECE

SAEVO INTRA SL 30

- Made of aluminium
- Intra-rotary connection system
- Weight: 50g
- Transmission: 1:1
- Drill Fixation: Friction Grip
- Drills: Rod type 2



External spray
adapter

WHEN CHOOSING YOUR HANDPIECES, BE AMAZED.

A high-performance academic kit, with completely redesigned handpieces!

ACADEMIC KIT INCLUDES:

- 1x Turbine
- 1x Contra-Angle
- 1x Micromotor
- 1x Straight Piece
- 1x Exclusive Bag



- Stronger and lighter parts
- Aluminium body
- Impact-resistant top
- Push-button system



TURBINE

Used to cut extremely rigid and resistant materials.

CONTRA-ANGLE

For its low-speed rotation. Used to remove softened dental structures, preserving the healthy tissue.

STRAIGHT PIECE

Ideal for finishing, polishing, and dental prosthetics.

MICROMOTOR



Functioning as a motor element for the counter-angle and the straight piece.

A COMPLETE HANDPIECE KIT - MEETING ALL THE DENTAL SURGEON'S NEEDS, AT OPTIMUM COST-BENEFIT. COMES IN A HANDY CARRY BAG.

ANVISA: TURBINE - 10069210065
ANVISA: MICROMOTOR - 10069210073
ANVISA: CONTRA-ANGLE - 10069210073
ANVISA: STRAIGHT PIECE - 10069210073

May 2023

REQUEST YOUR QUOTE
HANDPIECES



SAEVO is a member of the Alliage group. | Contact: sales@alliage-international.com

alliage