

DENTAL SET S 500 F

SAEVO

**TECHNOLOGY
ALIGNED WITH
YOUR TALENT**

The S500 brings you cutting-edge technology & design.



Illustrative image.

GET TO KNOW S500

- 1 Biarticulated headrest
- 2 Pneumatic patient arm support
- 3 Upholstery with 15 color options
- 4 Joystick control pedal
- 5 Anti-skid edging
- 6 No need for floor mounting
- 7 Coupled PAD/ Negatoscope kit
- 8 Allows coupling of up to 6 terminals
- 9 Electric Micromotor* Contra-Angle*
- 10 Manometer to visualize air pressure and power supply of the equipment
- 11 Support tray for instruments
- 12 Retractable rods with lock (version H)
- 13 Delivery unit with pneumatic arm
- 14 Flush system
- 15 180° folding vat with proximity sensor and cup-holder
- 16 Assistant arm with control pad
- 17 Dental light: Sirius 3x2 LEDs with proximity sensor
- 18 Multimedia support
- 19 Ergonomic stool*

[VIEW THE LIST OF STANDARD AND OPTIONAL ITEMS ON THE BACK](#)

*Optional.

DENTAL SET S 500 H



*Optional.

ELECTRIC MICROMOTOR*

With anti-noise technology, superior to pneumatic micromotors - ensuring total comfort for you and your patients. Experience high torque, precision, speed control, and productivity.



KEY FEATURES:

LOW NOISE

TORQUE

3,5 N.cm

ROTATION

- Multiplier contra-angle: 200000 rpm*
- Electric micromotor: up to 40000 rpm*

LIGHTING

25.000 lux

QUADRUPLE SPRAY

directed to the tip of the drill

*Optional.

SAEVO COLOURS

Upholstery that inspires your work routine

PVC OR VISCOELASTIC*



LEATHER*



*Optional. Illustrative images

STANDARD/OPTIONAL

See what comes with your **S 500**

N: NOT APPLICABLE | O: OPTIONAL | S: STANDARD ITEM

DENTAL SET

	F	H
Trendelenburg Position	S	S
Massage System	O	O
PVC Upholstery	S	S
Viscoelastic Upholstery	O	O
Leather Upholstery	O	O
Biarticulated Headrest	S	S
Armrests: 1 Fixed + 1 Folding and Pneumatic	S	S
Foot protector	S	S
Joystick Pedal (5 keys) - with Chip Blower Function	S	S
Foot Pedal - with Chip Blower Function (mobile)	O	O
Soft Comfort	O	O
Connection Box	O	O

DELIVERY UNIT

Pneumatic Lock Arm	S	S
Electrical Cables Kit	S	S
Upper Tray (Acrylic)	S	N
Lower Tray (Injected)	N	S
Kit PAD (Electronic Controls)	S	S
X-ray Viewer Kit	S	S
Flush Kit	S	S
Manometer	S	S
High Speed Terminal	S	S
Low Speed Terminal - with Spray	S	S
Injected Triple Syringe	S	S
Triple Syringe Heating Kit	O	O
Number of terminals in the delivery unit	UP TO 5 (3S & 2O)	UP TO 5 (3S & 2O)
Electric Micromotor Kit	O	O
Borden Terminal Kit - with Spray	O	O
Fibre Optic Terminal Kit	O	O
Curing Light Kit	O	O
Ultrasonic Scaler Kit	O	O
Bicarbonate Jet Kit	O	O

WATER UNIT

Cup Holder Kit - with Sensor	S	S
Cup Holder Kit - without Sensor	O	O
Electric Water Bowl Activator	S	S
Assistant Arm - with PAD control	S	S
Venturi Suction Kit	S	S
Vac Plus Suction Kit	S	S
External Suction Terminal Kit	O	O
Injected Triple Syringe Kit	O	O
Triple Syringe Heating Kit	O	O

DENTAL LIGHT

Sirius Sensor - 3 LEDs	O	O
Sirius Sensor - 5 LEDs	O	O
Sirius Sensor - 5 (3x2) LEDs	S	S

OTHERS

Multimedia Kit Support	O	O
------------------------	---	---

DENTAL SET S500 F AND S500 H CONFIGURATION:

Dental Chair: S500 • Anvisa: 10069210062

Water Unit: Syncrus G4A • Anvisa: 10069210063

Delivery Unit: Syncrus G8 FE / G8 HE / G3 / G3 H • Anvisa: 10069210075

Dental Light: Sirius G8 Sensor - 5 (3X2) LEDs • Anvisa: 10069210069

COD: 51000000216 • May 2023

REQUEST
YOUR QUOTE
S500