

THE

SAEVO SLIM

HAS ARRIVED

OUR NEW DIGITAL INTRAORAL SENSOR

Evolution that you see and your patient feels. Much more comfort in service with the image quality you already know.



+ THINNER
MORE COMFORT
FOR THE PATIENT

Saevo Slim offers images in **high definition** and connects your office with the digital world.

With it, you **optimize processes intraoral radiographic** capture and efficiently manages pictures of your patients.

BENEFITS

- Facilitates access among professionals
- Manipulates and save images
- Environmentally-friendly
- Provides accuracy in diagnosis
- Reduces film and processing costs

Size 2:
active area: 26x36mm

Size 1:
active area: 20x30mm

The Saevo Slim sensor is compatible with most equipment X-rays of the market, always **delivering high definition images**, with a great level of detail.

In addition, it has software that can be installed on an **unlimited number of computers** with tools for load, edit and share the images obtained.

Forget X-Ray films and developing liquids. save time and Contribute to the preservation of the environment!



FEATURES

UNLIMITED LICENSES

The Saevo Slim software can be installed on a unlimited number of computers. You can install it over a network, in a folder in the cloud, and access your radiographs from anywhere.

COMPATIBILITY

Saevo Slim Intraoral Sensor is compatible with existing intraoral X-ray equipment.

SHARE IMAGES AND DIAGNOSTICS

Save images in various formats like TIFF, BITMAP, PNG, DICOM, JPEG. Generate PDF reports and instantly share them with your clients and colleagues who will continue the treatment.

ORGANIZED IMAGES

Register your patients, capture images instantly and organize all images in seconds.



EASY CABLE REPLACEMENT

The Saevo Slim uses a USB 2.0 cable with 3m* and is easy to replace.

SHOCK RESISTANT

The internal protective silicon cover makes Saevo Slim shock resistant.

COMPLETE SOFTWARE

Wide range of tools such as rotation, zoom, brightness, contrast and gamma adjustment, positive/negative, color adjustment, embossing, text insertion, arrows, circles, area selection and cropping, angle and segment measurement, bar for reports/diagnosis, odontogram and comparison of full screen images

SAVE TIME

While it takes you about 5 minutes to capture and reveal an image through the analog system, Saevo Slim does it in about 5 seconds.

SAVE MONEY

Stop buying film developers and fixers. Instead of spending on the purchase of these products every month, invest in purchase the Saevo Slim. Generate results for you instead waste to the environment.

THE PERFECT PICTURE SHARPNESS

Capture high definition images in real time. With 25 line pairs (theoretical) and advanced image processing features, you have clear images for the best diagnosis.

THE SAEVO SLIM HAS THE CSI - TL

Technology Cesium Iodide (CsI) is the highest technology in High-Definition capture technology on the market. This technology, present today in tomographs, provides sharper images with reduced radiation radiation.

IMPRESS YOUR CUSTOMERS

Your clients will be impressed by the technology and richness of detail before and after the procedure.

ENVIRONMENTALLY FRIENDLY

The Saevo Slim eliminates the use of chemicals such as developers and fixers.

TECHNICAL INFORMATION

- Operational systems: Windows 7 Home, Start, Pro (x86 and x64); Windows 10 Home e Pro (x86 and x64)
- CPU: Intel Core i3 3.0 GHz or higher.
- RAM: 4GB
- HD: 500 GB or higher
- USB Connector: 2.0- Monitor 1024 X 768

SOLICITE ORÇAMENTO



Illustrative Photos. Saevo reserves the right to make changes to the design and technical specifications of the products upon authorization from ANVISA in order to continue improving them. Digital Sensor: Anvisa 10101139029. 51000000300 Version - Jan/2023. * variation of + or - 0.5m.