

## HARMONIZE YOUR TALENT

Saevo HOF is an dental set suitable for professionals exclusive in the area of cosmetic dentistry. It is ideal for those who follow market news in seeking to bring the best experience and comfort to their patients.

- LATERAL LIFT (COLUMN LIFT SYSTEM)
- SMALL BASE
- GREATER PROXIMITY TO THE OPERATIVE FIELD
- FOLDING LEG SUPPORT
- FRONTAL POSITIONING TO THE PATIENT
- VISCOELASTIC UPHOLSTERY
- BI-ARTICULATED HEADREST
- TRENDELENBURG POSITION
- JOYSTICK PEDAL



Illustrative picture. Check the list of options on the back.

### DENTAL CHAIR

- 1 Massage system \*
- 2 Leather Upholstery \*
- 3 Trendelenburg position.
- 4 Column Lift system
- 5 Viscoelastic Upholstery
- 6 Bi-articulated headrest
- 7 Folding armrest.
- 8 Folding leg support
- 9 Foot Protector

### 10 FOOT CONTROL

- Joystick, integrated into unit's base, with seat and headrest actuation

### 11 TRAY AND PANEL PAD

- **Tray:**
  - Auxiliary tray for instruments
  - Made of acrylic
- **Control pad:**
  - 3 programmable working positions, position back to zero, activation of the seat and backrest of the chair and activation of the reflector. \*

### DENTAL LIGHT \*

- 3 LED's / 5 LED's / 3x2 LED's \*
- Lower energy consumption and longer lifespan
- Proximity sensor
- Selection of luminous intensity
- Removable and autoclavable handles
- Removable front protector

\* Optional



# HOF

## SAEVO'S COLOURS

### Upholstery that inspires your work routine

#### VISCOELASTIC PVC



#### LEATHER\*



\*Optional. Illustrative picture

## CONFIGURATION

### Check out what comes with your HOF

O: OPTIONAL | S: STANDARD

| DENTAL CHAIR                            | SAEVO HOF |
|---|-----------|
| Column Lift System                      | S         |
| Trendelenburg Position                  | S         |
| Massage System                          | O         |
| Viscoelastic Upholstery                 | S         |
| Leather Upholstery                      | O         |
| Bi-articulated Headrest                 | S         |
| Chair Arm 1 Fixed + 1 Folding (Model V) | S         |
| Folding Leg Support                     | S         |
| Foot Protector                          | S         |
| Auxiliary tray kit**                    | O         |
| Soft Comfort                            | O         |
| Single Connection Box                   | O         |

| DENTAL LIGHT               | SAEVO HOF |
|----------------------------|-----------|
| Sirius Sensor 3 Leds       | O         |
| Sirius Sensor 5 Leds       | O         |
| Sirius Sensor 5 (3x2) Leds | O         |




| OTHERS            | SAEVO HOF |
|-------------------|-----------|
| Monitor support** | O         |

\*\*If you purchase the dental light

Illustrative photos. Colours may have changed in graphic reproduction. Saevo reserves the right to make changes to the design and technical specifications of the products upon authorization from ANVISA in order to continue improving them.  
CONFIGURATION: CHAIR: NEW VERSA - Anvisa 10069210062 | WATER UNIT: SYNCRUS G3 - Anvisa 10069210063  
DENTAL LIGHTS: Sirius G8 Sensor 3 LEDs / Sirius G8 Sensor 5 LEDs / Sirius G8 Sensor 5 (3X2) LEDs - Anvisa 10069210069

REQUEST YOUR QUOTE  
HOF

May 2023



SAEVO is a member of the Alliage group. | Contact: sales@alliage-international.com