

PORTABLE DENTAL X-RAYS AIR X-RAY



Illustrative image.



60 kV
2,5 mA

EVOLUTION IS IN YOUR HANDS

Mobility

Bringing efficiency and sophistication wherever you are.

Light and Compact

Light-weight and user-friendly, providing practicality that increases daily productivity.

Adaptable

Ready for use with an intraoral sensor, phosphorus plate and radiographic film.

SAEVO

AIR X-RAY



Illustrative Photos

- The Portable X-Ray ensures patient and operator safety, with mechanisms that reduce patient exposure.
- The Portable X-Ray minimizes radiation leakage via a radiopaque material collimator - avoiding unnecessary user exposure and directing the radiation to the exam target.
- Nominal voltage: 60kVp
Tube current: 2.5 mA
Focal spot: 0.7x0.7mm

TECHNICAL INFORMATION

- Ampoule Voltage:** 60kV
- Ampoule Current:** 2,5 mA
- Focus-skin Distance:** 200 mm
- Focus-receiver Distance:** 220 mm
- Maximum Energy Accumulated in 1h:** 150 mAs
- Weight:** 2,4 kg
- Maximum Output Power:** 150 W (60kV x 2,5mA)
- Battery:** Rechargeable Li-ion
- Generator:** Constant Power High Frequency Generator
- Target Angle:** 16°
- Maximum Operating Voltage (Intensity):** 60KVp
- Irradiation Time Selection Range:** 0.01 to 1 second (With steps of 0.01s)



ANVISA: 10069210087

May 2023

REQUEST
YOUR QUOTE
PORTABLE X-RAY